



Product Datasheet

Lactate Dehydrogenase Assay Kit (orb219863)

Description	Lactate Dehydrogenase Assay Kit
Reactivity	All
Range	30 µmol/L - 600 µmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Experiment Notes: Colorimetric
Sample Types	Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
CAS Number	9001-60-9
Expiration Date	12 months from date of receipt.



Representative
Standard
Curve



Positive
Control
reaction in 96-
well pla...