


Product Datasheet

NADP Assay Kit (orb219865)

Description	NADP Assay Kit
Reactivity	All
Range	0.1 µmol/L - 50 µmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Experiment Notes: Colorimetric
Sample Types	Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
CAS Number	53-59-8
Expiration Date	12 months from date of receipt.



Representative
Standard
Curve