


Product Datasheet

Hydrogen Peroxide Assay Kit (orb219868)

Description	Hydrogen Peroxide Assay Kit
Reactivity	All
Range	0.5 mmol/L - 100 mmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Experiment Notes: Colorimetric
Sample Types	Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
CAS Number	7722-84-1
Expiration Date	12 months from date of receipt.



Representative
Standard
Curve