


Product Datasheet

Trehalase Assay Kit (orb219871)

Description	Trehalase Assay Kit
Reactivity	All
Range	0.5 mmol/L - 5 mmol/L
Tested Applications	FA
Storage	Shipped and store at 4°C for 6 months.
Note	For research use only
Application notes	Experiment Notes: Colorimetric
Sample Types	Urine, Serum, Plasma, Tissue extracts, Cell lysate, Cell culture media, Other biological fluids
CAS Number	9025-52-9
Expiration Date	12 months from date of receipt.



Representative
Standard
Curve