


## Product Datasheet

### TRITC IF Detection System Kit (Mouse Anti-Goat) [Out of stock] (orb219883)

<b>Description</b>	TRITC IF Detection System Kit (Mouse Anti-Goat)
<b>Tested Applications</b>	ICC, IF
<b>Storage</b>	Store at 4°C for one year.
<b>Note</b>	For research use only
<b>Application notes</b>	Super Sensitive TRITC IF Detection System Kit is a non-biotin, one-step detection system that allows for the demonstration of antigens in paraffin-embedded tissue, cryostat sections, and cell preparations. This kit has been developed using a proprietary hyper labeling technology used to label IgG directly with more Tetramethylrhodamine (TRITC). One step polymer system provides increased sensitivity, time savings and detection simplicity. All the components contain with PBS, proteins, stabilizers and preservatives. This kit is suitable for single, double and triple immunofluorescence detect. The color is emerald green with a correct result. The TRITC is easily faded with light; all experiment process need keep away from light.
<b>Expiration Date</b>	12 months from date of receipt.



Immunofluorescent analysis staining in A...