

Product Datasheet

hsa-mir-21 Real-time RT-PCR Detection Kit [Out of stock] (orb219917)

Description	Detect mature miRNA hsa-mir-21 expression.
Tested Applications	PCR
Target	hsa-mir-21
Storage	Store at 4°C for short term and at -20°C for long term storage
Note	For research use only
Application notes	50 reactions for Real-Time PCR

Experiment Notes: The miRNAs Real-Time RT-PCR Detection Kit contains the whole reagent including the stem-loop RT primer and a pair of specific PCR primer set. It provides a sensitive and specific method of real-time PCR for quantificating of miRNAs from total RNA samples. It is Ideal for validating microRNA expression profiling data and high-throughput applications. Trait lies in even highly homogized miRNAs can be accurately dicriminated ensured by stemloop RT primer and high specific miRNAs primer set.

Sensitivity	10 copies of a standard mimics and as little as 25 pg of total RNA
Expiration Date	12 months from date of receipt.