

Product Datasheet

Anti-BCMA (human) (Vicky-1) Antibody (orb2202956)

Description	Anti-BCMA (human) (Vicky-1) Antibody
Species/Host	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FC, ICC, WB
Immunogen	Recombinant human BCMA extracellular domain (aa 2-54) (Prod. No. ALX-522-026).
Storage	Long Term Storage : -2°CShort Term Storage : -2°CHandling : Avoid freeze/thaw cycles.Stable for at least 1 year after receipt when stored as recommended.
Note	For research use only
Application notes	Application Info: ELISAFLOW CytometryImmunocytochemistryFunctional Application (blocks binding of BAFF to human BCMA)Optimal conditions must be determined individually for each application.
Isotype	IgG1
Clonality	Monoclonal
Purity	Protein G-affinity purified.
Uniprot ID	Q02223
Expiration Date	12 months from date of receipt.