

Product Datasheet

Goat Anti-TSPO / PBR (isoform) Antibody (orb2203710)

Description

Goat Anti-TSPO / PBR (isoform)

Species/Host	Goat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	TSPO
Preservatives	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Form/Appearance	Supplied at 0.5 mg/ml in Tris saline, 0.02% sodium azide, pH 7.3 with 0.5% bovine serum albumin. Aliquot and store at -20°C. Minimize freezing and thawing.
Storage	Aliquot and store at -2°C. Minimize freezing and thawing.
Note	For research use only
Application notes	Application Info: P-ELISA: Antibody detection limit dilution 1:128000.WB: Approx 17kDa band observed in Human Breast cancer lysates (calculated MW of 18.8kDa according to NP_000705.2). Recommended concentration: 0.1-0.3ug/ml.IHC: In paraffin embedded Human Colon shows strong cytoplasm staining of absorptive cells. Recommended concentration: 2-4ug/ml. Paraffin embedded Human Kidney. Recommended concentration 5 ug/mL.
Protein Sequence	C-RDNHGWRGGRRRLPE
Clonality	Polyclonal
MW	18.8
Uniprot ID	P30536
NCBI	NP_000705.2
Expiration Date	12 months from date of receipt.