

Product Datasheet

oar-mir-184 Real-Time RT-PCR Detection and U6 Calibration Kit [Out of stock] (orb220641)

Description

Detect mature miRNA oar-mir-

U6 snRNA.

Tested Applications

PCR

Target

oar-mir-184

Storage

Store at 4°C for short term and at -20°C for long term storage

Note

For research use only

Application notes

50 reactions for Real-Time PCR

Experiment Notes: The Real-Time RT-PCR detection and U6 calibration Kit contains the whole reagent including the stem-loop RT primer & specific PCR primer set for microRNA quantification and U6 primer set as housekeeping gene for calibration. It provides a sensitive and specific method of real-time PCR for relative quantification of miRNAs from total RNA samples. It is ideal for validating microRNA expression profiling data and high-throughput applications. Trait lies in even highly homogenized miRNAs can be accurately discriminated ensured by stemloop RT primer and high specific miRNAs primer set.

Sensitivity

10 copies of a standard mimics and as little as 25 pg of total RNA

Expiration Date

12 months from date of receipt.