



Product Datasheet

FAM113A antibody (orb2210)

Description	FAM113A antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human FAM113A (51-150/454 aa)
Target	FAM113A
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	52 kDa
Uniprot ID	Q9H1Q7
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Human sk...



Immunohistochemical staining of rat brai...