

Product Datasheet

IL16 Antibody (orb221297)

Description

IL16 Antibody

Species/Host

Rabbit

Reactivity

Human

Conjugation

Unconjugated

Tested

ELISA, ICC, IHC, WB

Applications

Immunogen

E.coli-derived human IL-16 recombinant protein (Position: S1212-S1332). Human IL-16 shares 86% amino acid (aa) sequence identity with mouse IL-16.

Form/Appearance

Lyophilized

Concentration

Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml.

Storage

Store at -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freeze-thaw cycles.

Note

For research use only

Application notes

WB: The detection limit for IL-16 is approximately 0.25ng/lane under reducing conditions. Tested Species: In-house tested species with positive results. By Heat: Boiling the paraffin sections in 10mM citrate buffer, pH6.0, for 20mins is required for the staining of formalin/paraffin sections. Other applications have not been tested. Optimal dilutions should be determined by end users. . Add 0.2ml of distilled water will yield a concentration of 500ug/ml.

Isotype

Rabbit IgG

Clonality

Polyclonal

MW

141752 MW

Uniprot ID



[Q14005](#)

Dilution Range

Immunohistochemistry (Paraffin-embedded Section), 0.5-1µg/ml, Human, By Heat Immunocytochemistry , 0.5-1µg/ml, Human, -Western blot, 0.1-0.5µg/ml, Human ELISA, 0.1-0.5µg/ml, Human

Expiration Date


12 months from date of receipt.



IHC(P)
analysis of
Human Tonsil
Tissue u...



WB analysis
of
Recombinant
Human IL-16
P...



WB analysis
of IL16 using
anti-IL16
anti...