

Product Datasheet

TNF alpha antibody (orb221324)

Description

TNF alpha antibody

Species/Host

Rabbit

Reactivity

Human, Mouse, Porcine, Rat

Conjugation

Unconjugated

Tested

ELISA, ICC, IF, IHC-P

Applications

Immunogen

KLH conjugated synthetic peptide derived from human TNF alpha. Please contact us for the exact immunogen sequence. The peptide is available as orb374840.

Target

TNF alpha

Form/Appearance

10 mM PBS, 0.02% sodium azide

Concentration

- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml

Storage

Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.

Note

For research use only

Isotype

IgG

Clonality

Polyclonal

Purity

Polyclonal antibodies are purified by peptide affinity chromatography

MW

17/26 kDa

Uniprot ID

[P06804](#), [P16599](#), [P01375](#)

NCBI

[000594](#), [000585](#)

Entrez



[7124](#)

Dilution Range


IF/ICC: 1:100, IHC-P: 1:100

Expiration Date


12 months from date of receipt.



Immunofluorescence analysis of pig liver...



IHC-P staining of pig liver tissue using...



IF image of pig liver vessel tissue usin...