



Product Datasheet

Snail antibody (orb221325)


Description	Snail antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human Snail. Please contact us for the exact immunogen sequence. The peptide is available as orb374841.
Target	Snail
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	29 kDa
Uniprot ID	Q95863, Q02085, Q6AY35
NCBI	21, 31, 005985, 005976
Entrez	6615
Dilution Range	IF/ICC: 1:100, IHC-P: 1:100, WB: 2 µg/ml
Expiration Date	12 months from date of receipt.



Western blot analysis of human breast ca...



IHC-P image of rat ihc-p tissue using an...



IHC-P image of pig ihc-p tissue using an...