



Product Datasheet

bFGF antibody (orb221347)


Description	bFGF antibody
Species/Host	Rabbit
Reactivity	Guinea pig, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ICC, IF, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human bFGF. Please contact us for the exact immunogen sequence. The peptide is available as orb393930.
Target	bFGF
Form/Appearance	10 mM PBS, 0.02% sodium azide
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml- 200 µg (in 400 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	30 kDa
Uniprot ID	P13109 , P15655 , P09038
NCBI	001997 , 002006
Entrez	2247
Dilution Range	IF/ICC: 1:100-600, IHC-P: 1:100-600, WB: 1:200-1000
Expiration Date	12 months from date of receipt.



IHC-P image of guinea pig lung tissue using...



WB analysis of mouse brain (lane 1) using...



IHC-P staining of rat colon tissue using...