

Product Datasheet

H5N6 hemagglutinin antibody (orb221432)

Description	H5N6 hemagglutinin antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from Influenza A virus H5N6 hemagglutinin (351-450/54 aa)
Target	H5N6 hemagglutinin
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	63 kDa
Expiration Date	12 months from date of receipt.