



Product Datasheet

Integrin alpha 5 antibody (orb221434)


Description	Integrin alpha 5 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human Integrin alpha-5 light chain (1021-1099/1099 aa)
Target	Integrin alpha 5
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	110 kDa
Uniprot ID	P08648
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of Mouse pl...



Immunohistochemical staining of human co...



Immunohistochemical staining of human lu...