

Product Datasheet

HOXD12 Antibody (orb2214672)

Description	HOXD12 Antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IF
Immunogen	Recombinant Human Homeobox protein Hox-D12 protein (30-169AA)
Preservatives	Liquid
Form/Appearance	Liquid
Storage	-20°C or -80°C
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	> 95%
Uniprot ID	P35452
NCBI	NP_067016
Expiration Date	12 months from date of receipt.