


Product Datasheet

MYO1H antibody (orb221548)

Description	MYO1H antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC, WB
Immunogen	KLH conjugated synthetic peptide derived from human MYO1H (501-600/958 aa)
Target	MYO1H
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	119 kDa
Uniprot ID	Q8N1T3
Dilution Range	WB=1:500-2000, ICC=1:100, IHC-P=1:100-500, IHC-F=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



Western blot: Lane 1:U87MG (Human) Cell;L...