


Product Datasheet

MX2 antibody (orb221565)

Description	MX2 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ICC
Immunogen	KLH conjugated synthetic peptide derived from human MX2 (621-715/715 aa)
Target	MX2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	79 kDa
Uniprot ID	P20591
Dilution Range	ICC=1:25
Expiration Date	12 months from date of receipt.



Immunohistochemical
staining of human
lu...