

Product Datasheet

NAP1/C12orf30 antibody (orb221601)

Description	NAP1/C12orf30 antibody
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human NAP1/C12orf30 (51-150/972 aa)
Target	NAP1/C12orf30
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	112 kDa
Uniprot ID	Q14CX7
Expiration Date	12 months from date of receipt.