

Product Datasheet

NDUFB1 antibody (orb221669)

Description	NDUFB1 antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human NDUFB1 (21-58/58 aa)
Target	NDUFB1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	7 kDa
Uniprot ID	075438
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.