



Product Datasheet

NOTCH4 antibody (orb221726)


Description	NOTCH4 antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthetic peptide derived from human NOTCH4 (1151-1250/2003 aa)
Target	NOTCH4
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	207 kDa
Uniprot ID	Q99466
Dilution Range	IHC-P=1:100-500
Expiration Date	12 months from date of receipt.



IHC-P
analysis of
Mouse brain
using NOTC...



IHC-P
analysis of
Rat liver
using
NOTCH4...



IHC-P
analysis of
human
myocardium
using...