

Product Datasheet

NKF3 kinase family member antibody (orb221756)

Description	NKF3 kinase family member antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	KLH conjugated synthetic peptide derived from human NKF3 kinase family member (661-760/1746 aa)
Target	NKF3 kinase family member
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	193 kDa
Uniprot ID	Q9H792
Expiration Date	12 months from date of receipt.