

Product Datasheet

GRIN2A (phospho-Ser1232) antibody (orb221779)

Description	GRIN2A (phospho-Ser1232) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human NMDAR2A around the phosphorylation site of Ser1232 (MR(p-S)PF)
Target	phospho-NMDAR2A (Ser1232)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	163 kDa
Uniprot ID	Q12879
Expiration Date	12 months from date of receipt.