


Product Datasheet

nNOS (phospho-Ser847) antibody (orb221792)

Description	nNOS (phospho-Ser847) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from mouse nNOS around the phosphorylation site of Ser847 (FN(p-S)VS)
Target	phospho-nNOS (Ser847)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	160 kDa
Uniprot ID	Q9Z0J4
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of human lu...