


Product Datasheet


Neural retinal specific leucine zipper antibody (orb221839)

Description

Neural retinal specific leucine zipper



Immunohistochemical staining of Rat eye ...



Western blot analysis of human Raji cell...

Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	KLH conjugated synthetic peptide derived from human Neural retinal specific leucine zipper/NRL (151-237/237 aa)
Target	Neural retinal specific leucine zipper/NRL
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid: 0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	affinity purified by Protein A
MW	26kDa
Uniprot ID	P54845
Dilution Range	WB=1:500-2000
Expiration Date	12 months from date of receipt.