

Product Datasheet

PLB (phospho-Thr17) antibody (orb221859)

Description	PLB (phospho-Thr17) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human PLB around the phosphorylation site of Thr17 (AS(p-T)IE)
Target	phospho-PLB (Thr17)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	6 kDa
Uniprot ID	P26678
Expiration Date	12 months from date of receipt.