



Product Datasheet

Paxillin (phospho-Ser83) antibody (orb221861)

Description	Paxillin (phospho-Ser83) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated Synthesised phosphopeptide derived from human Paxillin around the phosphorylation site of Ser83 (PQ(p-S)SS)
Target	Phospho-Paxillin (Ser83)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	68 kDa
Uniprot ID	P49023
Expiration Date	12 months from date of receipt.



Western blot analysis of mouse lung tiss...



Immunohistochemical staining of human la...