

Product Datasheet

RUNX1 (phospho-Ser276) antibody (orb221879)

Description	RUNX1 (phospho-Ser276) antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	KLH conjugated synthesised phosphopeptide derived from human RUNX1 around the phosphorylation site of Ser276 (PI(p-S)PG)
Target	phospho-RUNX1 (Ser276)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	49 kDa
Uniprot ID	Q01196
Expiration Date	12 months from date of receipt.