

Product Datasheet

Human Fab Of Normal IgG (orb22193)

Description Chromatographically purified Fab fragment of human IgG. Purified human Fab may be used as a referenc...

Tested Applications DOT, ELISA, WB

Storage Storage: After arrival store at +4°C. Prolonged storage may be at -20°C. Allow to stand at ambient temperature for 5-10 minutes to reach equilibrium. Reconstituted solution may be stored frozen. Working dilutions may be prepared by adding the required amount of phosphate buffered saline (PBS), should not be refrozen and preferably used within 24 hours. Repeated thawing and freezing should be avoided. If a slight precipitation occurs upon storage, this should be removed by centrifugation and will not affect the performance of the product

Note For research use only

Hazard Information This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals. This datasheet is as accurate as reasonably achievable, but Nordic-MUBio accepts no liability for any inaccuracies or omissions in this information.

Expiration Date 6 months from date of receipt.