

Product Datasheet

GPR55 antibody (FITC) (orb222085)

Description	GPR55 antibody (FITC)
Species/Host	Rabbit
Reactivity	Canine, Equine, Gallus, Human, Mouse, Rat
Conjugation	FITC
Tested Applications	ICC, IF
Immunogen	Synthetic peptide (please contact us for details)
Target	GPR55 (FITC)
Form/Appearance	Liquid: PBS, 2 mM Sodium Azide, 30% glycerol
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	37 kDa
Uniprot ID	Q9Y2T6
Entrez	9290
Dilution Range	IF/ICC:1:50-200 (based on 0.5 mg/ml)
Expiration Date	12 months from date of receipt.