



Product Datasheet

Integrin alpha 5 antibody (FITC) (orb222105)


Description	Integrin alpha 5 antibody (FITC)
Species/Host	Rabbit
Reactivity	Human, Mouse, Porcine, Rat
Conjugation	FITC
Tested Applications	ICC, IF
Immunogen	Synthetic peptide (please contact us for details)
Target	Integrin alpha 5 (FITC)
Form/Appearance	Liquid: PBS, 2 mM Sodium Azide, 30% glycerol
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	110 kDa
Uniprot ID	P08648, P11688
Entrez	3678
Dilution Range	IF/ICC:1:50-500 (based on 0.5 mg/ml)
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of mouse uterus.



Immunofluorescence image of pig large intestine.



IF analysis of pig large intestine tissue.