

Product Datasheet

MMP9 antibody (FITC) (orb222128)

Description	MMP9 antibody (FITC)
Species/Host	Rabbit
Reactivity	Bovine, Canine, Equine, Guinea pig, Human, Mouse, Porcine, Rat
Conjugation	FITC
Tested Applications	ICC, IF
Immunogen	Synthetic peptide (please contact us for details)
Target	MMP9 (FITC)
Form/Appearance	Liquid: PBS, 2 mM Sodium Azide, 30% glycerol
Concentration	- 100 µg (in 200 µl): 0.5 mg/ml
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Polyclonal antibodies are purified by peptide affinity chromatography
MW	78 kDa
Uniprot ID	P50282, P14780, P41245
Entrez	4318
Dilution Range	IF/ICC:1:50-200 (based on 0.5 mg/ml)
Expiration Date	12 months from date of receipt.