

Product Datasheet

Mast Cell Chymase antibody (FITC) [Out of stock] (orb222129)

Description	Mast Cell Chymase antibody (FITC)
Species/Host	Rabbit
Conjugation	FITC
Immunogen	KLH conjugated synthetic peptide derived from human CMA1
Target	CMA1
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P21844
Expiration Date	12 months from date of receipt.