

Product Datasheet

TMEM14A Antibody (orb2224592)

Description	TMEM14A Antibody
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	Recombinant Human Transmembrane protein 14A protein (45-78AA)
Preservatives	Liquid
Form/Appearance	Liquid
Storage	-20°C or -80°C
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	Q9Y6G1
NCBI	NP_054770
Expiration Date	12 months from date of receipt.