



Product Datasheet

CD64 antibody (Biotin) (orb222723)

Description	Mouse monoclonal antibody against CD64 conjugated to Biotin
Reactivity	Human, Primate
Conjugation	Biotin
Tested Applications	FC, ICC, IHC-Fr, IHC-P, IP, WB
Immunogen	Rheumatoid synovial fluid cells and fibronectin purified human monocytes
Target	CD64
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 4 µg/ml.
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Purity	Purified antibody is conjugated with biotin LC-NHS ester under optimum conditions and unconjugated antibody and free biotin are removed by size-exclusion chromatography.
Uniprot ID	P12314
Entrez	2209
Dilution Range	Flow cytometry: Recommended dilution: 4 µg/ml.
Expiration Date	12 months from date of receipt.



Separation of human monocytes (red-fille...



Flow cytometry surface staining pattern ...