


Product Datasheet

Helios antibody (PE) (orb222727)

Description

Hamster monoclonal antibody against Helios conjugated



Human peripheral blood lymphocytes were ...

Reactivity

Human, Mouse

Conjugation

PE

Tested Applications

FC

Immunogen

Peptide corresponding to the amino acids 51-107 of Helios

Target

Helios

Preservatives

Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide

Storage

Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.

Note

For research use only

Application notes

Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests. Intracellular staining.

Isotype

Armenian Hamster IgG

Clonality

Monoclonal

Purity

Purified antibody is conjugated with R-phycoerythrin (PE) under optimum conditions. Unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.

Uniprot ID
[Q9UKS7](https://www.uniprot.org/entry/Q9UKS7)
Entrez
[22807](https://www.ncbi.nlm.nih.gov/entrez/query.fcgi?db=Gene&term=22807)
Dilution Range

Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests. Intracellular staining.

Expiration Date

12 months from date of receipt.