


Product Datasheet

CD253 antibody (orb222728)

Description	Mouse monoclonal antibody to CD253
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FA, FC
Immunogen	Recombinant soluble fragment (aa 95-281) of human TRAIL.
Target	CD253
Preservatives	Phosphate buffered saline (PBS), pH 7.4
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Low Endotoxin, Azide free format. Product does not contain preservatives and is intended primarily for functional applications. Preservative containing buffers are available for IHC applications
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P50591
Entrez	8743
Dilution Range	Functional application: The antibody 2E5 has high neutralizing activity for human TRAIL in biological assays. Flow cytometry: Recommended dilution: 1-10 µg/ml.
Expiration Date	12 months from date of receipt.



Functional application example: Apoptos...