



Product Datasheet

CD32 antibody (orb222729)


Description	Mouse monoclonal antibody to CD32
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, MC
Immunogen	purified glycosylated recombinant human FcγRIIa2
Target	CD32
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 2-6 µg/ml.
Isotype	Mouse IgG1 kappa
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P12318
Entrez	2212
Dilution Range	Flow cytometry: Recommended dilution: 2-6 µg/ml.
Expiration Date	12 months from date of receipt.



Mass cytometry (surface staining) of hum...



Separation of human CD32 positive lympho...



Flow cytometry surface staining pattern ...