

Product Datasheet

LRBA Antibody - HRP Conjugated (orb2227343)

Description	LRBA Antibody - HRP Conjugated
Species/Host	Rabbit
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Human Lipopolysaccharide-responsive and beige-like anchor protein (2-252AA)
Preservatives	Liquid
Form/Appearance	Liquid
Storage	-20°C or -80°C
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Caprylic Acid Ammonium Sulfate Precipitation purified
Uniprot ID	P50851
NCBI	NP_006717.2, NP_001186211.2
Expiration Date	12 months from date of receipt.