



Product Datasheet

Sialoadhesin antibody (orb222739)

Description	Mouse monoclonal antibody to Sialoadhesin
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-Fr, IP, WB
Immunogen	human rhinovirus 14-infected monocyte-derived dendritic cells
Target	CD169
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-4 µg/ml
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	Q9BZZ2
Entrez	6614
Dilution Range	Flow cytometry: Recommended dilution: 1-4 µg/ml
Expiration Date	12 months from date of receipt.



Separation of human monocytes (red-fille...



Flow cytometry surface staining pattern ...