

## Product Datasheet

### Mouse Mab Human IgG3 (FITC) antibody (orb22274)

**Description**

Fluorescein isothiocyanate-conjugated

**Species/Host**

Mouse

**Conjugation**

FITC

**Preservatives**

Purified monoclonal mouse IgG1 kappa conjugated with FITC, lyophilized from a solution in phosphate buffered saline (pH7.2). No preservative added, as it may interfere with the antibody activity. No foreign protein added.

**Form/Appearance**

Purified monoclonal mouse IgG1 kappa conjugated with FITC, lyophilized from a solution in phosphate buffered saline (pH7.2). No preservative added, as it may interfere with the antibody activity. No foreign protein added.

**Storage**

Storage: The lyophilized product is shipped at ambient temperature and may be stored at +4°C; prolonged storage at or below -20°C. Dilutions may be prepared by adding phosphate buffered saline (PBS, pH 7.2). Avoid repeated thawing and freezing. If a slight precipitation occurs upon storage, this should be removed by centrifugation and will not affect the performance of the product. Diluted solutions should be stored at +4°C, not refrozen, and preferably used the same day

**Note**

For research use only

**Isotype**

IgG1

**Clonality**

Monoclonal

**Source**

Highly purified monoclonal IgG3m(u) isolated from human serum.

**Hazard Information**

This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals. This datasheet is as accurate as reasonably achievable, but Nordic-MUBio accepts no liability for any inaccuracies or omissions in this information.