

Product Datasheet

Sialoadhesin antibody (PE) (orb222740)

Description

Mouse monoclonal antibody against Sialoadhesin

Reactivity

Human

Conjugation

PE

Tested

FC

Applications

Immunogen

human rhinovirus 14-infected monocyte-derived dendritic cells

Target

CD169

Preservatives

Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide

Storage

Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.

Note

For research use only

Application notes

Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 10⁶ cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests.

Isotype

Mouse IgG1

Clonality

Monoclonal

Purity

Purified antibody is conjugated with R-phycoerythrin (PE) under optimum conditions. Unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.

Uniprot ID

[Q9BZZ2](#)

Entrez


[6614](#)

Dilution Range

Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 10⁶ cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests.

Expiration Date

12 months from date of receipt.



Flow
cytometry
analysis
(surface
stainin...