

Product Datasheet

CD129 antibody (orb222742)

Description	Mouse monoclonal antibody to CD129
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC
Immunogen	Human CD129-transfected cell line
Target	CD129
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: It is recommended to use bright fluorochromes or signal multiplying detection.
Isotype	Mouse IgG2b kappa
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	Q01113
Entrez	3581
Dilution Range	Flow cytometry: It is recommended to use bright fluorochromes or signal multiplying detection.
Expiration Date	12 months from date of receipt.