

Product Datasheet

CD129 antibody (orb222743)

Description	Mouse monoclonal antibody to CD129
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FA, FC
Immunogen	Human CD129-transfected cell line
Target	CD129
Preservatives	Phosphate buffered saline (PBS), pH 7.4
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Low Endotoxin, Azide free format. Product does not contain preservatives and is intended primarily for functional applications. Preservative containing buffers are available for IHC applications
Isotype	Mouse IgG2b kappa
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	Q01113
Entrez	3581
Dilution Range	Functional application: Blocking. Flow cytometry: It is recommended to use bright fluorochromes or signal multiplying detection.
Expiration Date	12 months from date of receipt.