



Product Datasheet

CTLA4 antibody (orb222749)


Description	Mouse monoclonal antibody to CTLA4
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, ICC, IHC-Fr, IP
Immunogen	Human CD152-IgG heavy chain fusion protein
Target	CD152
Preservatives	Phosphate buffered saline (PBS), pH 7.4
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Low Endotoxin, Azide free format. Product does not contain preservatives and is intended primarily for functional applications. Preservative containing buffers are available for IHC applications
Isotype	Mouse IgG2a
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	P16410
Entrez	1493
Dilution Range	Flow cytometry: Recommended dilution: 1-4 µg/ml, Intracellular staining.
Expiration Date	12 months from date of receipt.



Flow
cytometry
multicolor
surface
staini...



Separation of
human
CD152
positive CD3
p...



Flow
cytometry
surface
staining
pattern ...