


Product Datasheet

CD255 antibody (orb222754)

Description	Mouse monoclonal antibody to CD255
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC
Immunogen	human CD255-transfected 2PK-3 cells
Target	CD255
Preservatives	Tris buffered saline (TBS), pH 8.0, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Isotype	Mouse IgG3 kappa
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	O43508
Entrez	8742
Dilution Range	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Expiration Date	12 months from date of receipt.



Separation of
CD255
transfected
I5178Y c...