


Product Datasheet

CD255 antibody (orb222755)

Description	Mouse monoclonal antibody to CD255
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FA, FC, IHC
Immunogen	human CD255-transfected 2PK-3 cells
Target	CD255
Preservatives	Tris buffered saline (TBS), pH 8.0
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Low Endotoxin, Azide free format. Product does not contain preservatives and is intended primarily for functional applications. Preservative containing buffers are available for IHC applications
Isotype	Mouse IgG3 kappa
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	043508
Entrez	8742
Dilution Range	Functional application: Neutralization. Flow cytometry: Recommended dilution: 1-4 µg/ml
Expiration Date	12 months from date of receipt.



Separation of
CD255
transfected
I5178Y c...