


Product Datasheet

CD79b antibody (APC) (orb222757)

Description

Mouse monoclonal antibody against CD79b conjugated


Reactivity

Human

Conjugation

APC

Tested Applications

FC

Immunogen

Fraction of Ig-associated molecules isolated from Ramos B cells

Target

CD79b

Preservatives

Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide

Storage

Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.

Note

For research use only

Application notes

Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests.

Isotype

Mouse IgG1 kappa

Clonality

Monoclonal

Purity

Purified antibody is conjugated with activated allophycocyanin (APC) under optimum conditions and unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.

Uniprot ID

[P40259](#)

Entrez

[974](#)

Dilution Range

Flow cytometry: The reagent is designed for analysis of human blood cells using 10 µl reagent / 100 µl of whole blood or 106 cells in a suspension. The content of a vial (1 ml) is sufficient for 100 tests.

Expiration Date

12 months from date of receipt.

Flow
cytometry
analysis
(surface
stainin...