

Product Datasheet

Ovotransferrin Antibody - HRP Conjugated (orb2227584)

Description	Ovotransferrin Antibody - HRP Conjugated
Species/Host	Rabbit
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Chicken Ovotransferrin (Native Protein)
Preservatives	Liquid
Form/Appearance	Liquid
Storage	-20°C or -80°C
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	> 95%
Uniprot ID	P02789
Expiration Date	12 months from date of receipt.